

Contaminant Detection in Soybean Using Allele-specific Digital PCR on the TaqMan® OpenArray® Platform



Marion Webster¹, Wing Cheung², Sunali Patel³, Kevin Munnely³, Elena Grigorenko³

¹Life Technologies, 850 Lincoln Center Drive, Foster City, CA 94404, ²DNA LandMarks Inc., St-jean-sur-richelieu, J3B 6X3 Quebec, Canada, ³Life Technologies, 12 Gill Street, Woburn, MA 01801

Abstract

Contamination of seed lots by unexpected varieties is a problem in quality assurance of certified seeds. Identification and elimination of contaminated seed lots are essential to maintain the quality of seed stocks. Cost efficient and high-throughput methods for identifying low contamination level are needed.

We are presenting a feasibility study for contaminant detection in soybean using digital PCR on the TaqMan® OpenArray® platform. Digital PCR is a sensitive approach to detect rare alleles and allows absolute quantification without using endogenous controls. Combining digital PCR with the TaqMan® OpenArray® system enables high-throughput contaminant screening.

Six SNPs were selected to distinguish soybean variety A from the contaminant variety B. Allele specific Custom TaqMan® Assays were designed for the OpenArray® digital PCR. A spike-in experiment was set up mixing soybean variety A DNA with contaminant variety B DNA at a 10,000:1 ratio to mimic 0.01% contamination in the seed lot. Digital PCR was performed on the OpenArray® PCR instrument using real-time signal acquisition followed by data analysis for quantification of detected copies/μl of samples.

Specificity of the assays in discriminating contaminant in a background of high concentration of variety A DNA was tested. All contaminant-specific assays produced no false-positive results in variety A DNA background and achieved the targeted sensitivity of detecting 0.01% contamination.

This methodology combining digital PCR with the high-throughput OpenArray® platform is applicable to any type of rare allele detection like contaminant detection or identifying genetically modified plant seeds in a pool of wild-type seeds.

The Experiment: Digital PCR

Custom TaqMan® assays for detection of variety A and contaminant variety B were designed and validated using gDNA. Contaminant variety B was spiked into variety A at a 1:5,000 and 1:10,000 ratio. We used TaqMan® OpenArray® Digital PCR Plates for detection of contaminant variety B.




Figure 2. Digital PCR workflow. Samples are partitioned into individual through-holes at a concentration of 1 copy/μl or lower.

The Result: Sensitive Detection

The experiment was run using TaqMan® OpenArray® Digital PCR Plates. The OpenArray Digital PCR software was used for data analysis. The amplification heat map (below) shows through-holes with detected (red) or undetected (black) target amplification. The DNA copy number was calculated using Poisson analysis by counting the number of wells with detected amplification.




Figure 4. Heat Map showing a non-specific assay (left) amplifying both variety A and variety B alleles and a specific assay (right) amplifying its target variety B only. Sensitivity level: 1:10,000.

The Problem: Contamination

This study aims at detecting contaminant variety B in a population of variety A in the soybean DNA. Targeted sensitivity level of detection is 1:10,000. We mimicked this ratio by spiking contaminant variety B gDNA into variety A gDNA background.




Figure 1. Simulating seed contamination with a spike-in experiment, adding small amounts of contaminant DNA

The Platform: The OpenArray®

The OpenArray® platform offers high-throughput flexibility for different genomics application including digital PCR.




Figure 3. OpenArray® instrument and detail of the OpenArray® plate showing one of 3,072 through-holes

Summary

- We designed and validated 6 allele specific TaqMan® assays for SNPs that distinguish our strains of interest.
- Using digital PCR, we achieved the desired detection sensitivity of 1:10,000 contaminant variety B in variety A soybean seed DNA background.
- The digital PCR can be used for variety of AgBio applications including GMO testing.

References:

- Vogelstein B., Kinzler, KW. *Digital PCR*, PNAS, 1999, v.96, 9236-9241
- Morrison T, et al. *Nanoliter high-throughput quantitative PCR*. Nucleic Acids Res., 2006, v.34, e123